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
GB INSTRUCTIONS SUPER KEGAL MUSCLETRAINER

The unique design of the Super Kegel makes it easy for anyone to use. It's important however to position the exerciser in a way; that's comfortable and creates the leverage needed to close the wings to their fullest capacity.

Before you begin, be sure you are wearing comfortable clothes such as sweats, leotards etc., or for men supportive underwear sure that the tension device is set on the lowest tension setting. First, while holding the front saddle of the Super Kegel, place it between your upper thighs an inch or two below the pelvic illustrated. Do not position it directly against the pelvic floor as it will not allow the proper leverage to close the wings together it may pinch. When you have the Super Kegel in this position, while still holding the front saddle, pull it forward until you feel the wings snug against the back of your thighs. This positioning should afford the greatest amount of leverage for your repetition maximize development. Once you have the Super Kegel properly positioned, refer to the exercise chart for sets and repetitions. Be aware of your posture when you do the exercises. Do not lean to one side or tend a knee in towards the other while trying to close the wings. Keep your pelvis Straight and square with your knees slightly bent and pointing Straight ahead throughout your repetitions. Once you close the wings together and complete your sets and repetitions, you can reset the tension device to create even greater resistance for increased development. Touching the wings together may be difficult at first and you may have to work up to it. This should be an indication that you need to develop more strength in this entire area. You can monitor your improvement by using width of your finger(s) to measure how far you can close the wings as the program progresses, it's important to close the wings together because even though the initial pressure by the gluteal muscle brings the PO muscle into a flexed condition, it's the continued applied resistance to close the wings that creates the residual development to the PO muscle and shortens development time.

EU: MediaShop GmbH | Schneiderstraße 1, Top 11 2620 Neunkirchen | Austria | DE | AT | CH: OvSOO 376 36 01 office@mediashop-group.com | www.mediashop.tv

22203 | 102022 | Made in China

Forgalmazza:
Telemarketing International Kft.
H-9028 Győr, Fehérvári út 75.
Tel.: + 36-96/961-000
Fax: + 36-96/616-662
www.wsteshop.hu