# VibroSculpt

Gebrauchsanleitung .....	<b>02</b>
Instructions for Use .....	<b>05</b>
Mode d'emploi .....	<b>07</b>
Manuale di utilizzazione .....	<b>10</b>
Gebruiksaanwijzing .....	<b>13</b>
Használati utasítás .....	<b>16</b>
Návod k použítí .....	<b>19</b>
Návod na použitie .....	<b>22</b>
Instrucțiuni de utilizare .....	<b>24</b>

AS SEEN ON  
**TV**  
ORIGINAL

## **SICHERHEITS- INFORMATION**

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Hergestellt in China

## SAFETY INFORMATION

- Keep the device and its packaging away from children.
- Only operate using the supplied AC adapter (BYX-1202000E).
- Do not use if there is any damage evident upon the device or accessories.
- Keep the device, power supply and cables away from water, steam and other liquids.
- Never use in a bathroom, sauna, etc.
- Never attempt to recover the device if it falls into water. Disconnect the plug immediately.
- Never attempt to repair the device yourself. Incorrect repairs can pose considerable danger to the user.
- Switch off and unplug the device after each use and before each cleaning.
- This device may be used by children from 8 years of age, or by persons with reduced physical, sensory or mental capabilities or lack of experience and / or knowledge, if they have been supervised or instructed in the safe use of the device and have understood the hazards arising therefrom. Children are not allowed to play with the machine. Cleaning and user maintenance must not be performed by children unless they are over 8 years old and supervised. Children younger than 8 years of age are to be kept away from the machine and the connecting cable. Extreme caution is advised when using a device near children and when leaving it unattended.
- Do not use the product if you suffer from any ailments or if you are pregnant.
- Do not massage body areas that are swollen, inflamed, or have a rash.
- If in doubt, consult your doctor before using the product.
- Read the manual carefully.
- Use the product only for its intended purpose as described in the manual.
- This product is not intended to act as an alternative to professional medical advice, diagnosis, or treatment.

- You use this product at your own risk. The manufacturer and the distributor assume no liability for the use of this product.

### **CORRECT USE**

This device is solely intended for cosmetic massages and home use; it is not intended for medical or professional use. It should not be used as an alternative to medical treatment. Please seek medical advice if you are uncertain as to whether the massager is suitable for you. Do not use the device directly on varicose veins or following vein treatment or surgery. In such cases, we recommend that you always seek medical advice before using the device. You can use the treatment as often as necessary. To maintain the desired result on your legs, continue to use at least once a month.

### **PRECAUTIONARY MEASURES**

- This is not a medical device, and is exclusively intended for massages. No medical claims are implied or guaranteed by the use of this product.
- Ensure that the unit is turned off before connecting the mains adapter.
- It is recommended that you first start the massager at low speed/intensity and then gradually increase the speed/intensity as desired.
- Do not massage a single area for any longer than 10 minutes.
- Do not leave the product in places where children can access it.
- Do not use on varicose veins or damaged skin.
- Be careful when applying to birthmarks. If necessary, consult your doctor.
- Do not use on visible or fresh scars.
- Avoid running directly over the navel. Desist from using the device if the stomach is hard and hot, and consult a doctor.
- Possible redness of the treated skin is normal and only temporary.
- If the VibroSculpt seems to be faulty, do not

attempt to open it up or repair it.

- DO NOT USE WHEN PREGNANT OR IF STILL BREASTFEEDING. As with any massage product, people with special medical conditions are advised to seek medical advice before using.
- KEEP AWAY FROM LONG HAIR, JEWELLERY AND LOOSE CLOTHING.
- DO NOT USE ON THE HEAD, FACE, NECK, SPINE OR CHEST AREAS.

### **USING THE DEVICE:**

1. The VibroSculpt Cosmetic Massager is ONLY designed for use with the special massage pad that should be applied to the skin during every treatment.
2. Ensure that the switch is in the OFF position.
3. Plug the power adapter into the outlet. Make sure that the voltage indicated on the adapter corresponds to the voltage of the outlet.
4. Select the pad for the area of the body to be treated:  
BLUE for a light, smooth and gentle massage  
WHITE for a relaxing and dry massage  
ORANGE to reduce the appearance of cellulite
5. Ensure that the device is turned off before you exchange massage cloths. Then you can fit the desired pad.
6. Customize your massage by increasing or decreasing the speed.

### **USING THE SPECIAL MASSAGE PADS**

- The blue extra soft cloth is perfect for use on your whole body: it smooths, firms and tones.
  - The white cloth is ideal for relieving tension and for dry massage of the back, shoulders and legs.
  - The orange ultra-peeling cloth stimulates the skin to reduce the appearance of cellulite.
- #### **MASSAGING THE ABDOMEN AND HIPS**
- Lie on a bed or stand in a comfortable position. Keep your back straight and your knees bent. On a low-intensity setting, move the device away from the navel and across the abdomen. Then move sideways, working on the hip area on either side for a few minutes.

- The device can also be started from the centre of the abdomen towards the sides in straight lines. Repeat this motion several times, selecting a different starting point each time.
- This massage should only last for a short duration; never longer than 5-10 minutes.

### **MASSAGING THE ABDOMEN AND HIPS**

- Lie on a bed or stand in a comfortable position. Keep your back straight and your knees bent. On a low-intensity setting, move the device away from the navel and across the abdomen. Then move sideways, working on the hip area on either side for a few minutes.
- The device can also be started from the centre of the abdomen towards the sides in straight lines. Repeat this motion several times, selecting a different starting point each time.
- This massage should only last for a short duration; never longer than 5-10 minutes.

### **MASSAGING THE LEGS**

- To massage your legs, work from the bottom to the top. While sitting or lying down, allow the device to work on the insides and outsides of the knees for a few minutes. Continue on both sides of the thighs to the groin and hips.
- Let the device run for 5-10 minutes, always moving slowly upwards.
- Do not apply pressure to the device. All motions should be slow. Avoid fast, inaccurate movements.

### **MASSAGING THE BUTTOCKS AND HIPS**

- Move the device upward from the lower buttocks until you reach the hips.
- The best position for this is lying on the stomach, ensuring you have enough space for the application.

### **MASSING THE BACK**

- Sit on a bed or sofa, or stand in a comfortable position. Keep your back straight. Massage from the centre outwards. Put the device on a

low setting, and move it outwards in a sideways direction. Run the appliance for a few minutes on each side of the back.

- The device can also be started from the centre of the back towards the sides in straight lines. Repeat this motion several times, selecting a different starting point each time.
- This massage should be short, just 5-10 minutes.
- **DO NOT APPLY TO THE NECK AREA AND KEEP AWAY FROM LONG HAIR, JEWELLERY AND LOOSE CLOTHING. DO NOT APPLY TO THE HEAD, FACE, NECK, SPINE OR CHEST AREAS.**

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**EN:** Do not dispose this product at the end of its service life through normal domestic waste. It has to be disposed through a collection agency that recycles electric and electronic devices. This is pointed out by the symbol on the product, in the instruction manual, and on the packaging. Please inquire about local collection agencies that may be operated by your distributor or your municipal administration. By recycling, making use of the

materials, or other forms of recycling old units you are making an important contribution to protecting our environment.

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# Meal Plan

	BREAKFAST	SNACK	LUNCH	SNACK	DINNER
<b>MONDAY</b>	1 apple 4 rice cakes	Buttermilk (0.1% fat)	Rice with Chicken	1 tangerine 5 nuts	Tuna salad
<b>TUESDAY</b>	1 boiled egg 1 whole grain bread	2 apples	Steak with beans	1 banana	Vegetable soup
<b>WEDNESDAY</b>	Coffee or unsweetened tea 4 graham crackers	Buttermilk (0.1% fat) with berries	Grilled salmon with lemon rice	Yogurt with cucumber	Sautéed mushrooms with egg
<b>THURSDAY</b>	1 yogurt 1 orange	Raw vegetables	Vegetable casserole	Homemade smoothie	Tuna omelette
<b>FRIDAY</b>	Whole grain bread, nonfat farmers cheese	4 rice cakes with low fat topping	Garden salad with grilled chicken breast	Melon with prosciutto	Tomatoes with mozzarella cheese
<b>SATURDAY</b>	Farmers cheese with berries	Celery with yoghurt dip	Pikeperch filet with vegetables	1 whole grain bread with avocado	Salad greens with fried egg
<b>SUNDAY</b>	Pancakes with fresh fruit	1 tangerine 5 nuts	Stuffed peppers with tomato sauce	1 cucumber 5 nuts	Pumpkin-carrot soup

Drink 2 liters of water daily. Season your food with fresh herbs. Chili, ginger and cinnamon help you lose weight. Fry meats and vegetables without oil in a nonstick pan. Marinate salad with vinegar and oil.

Repeat this meal plan for two to three weeks. Vary the meals by alternating their order. Fruits and vegetables are interchangeable: For example, you can eat an apple instead of a tangerine.

Advice: Please talk to your doctor before starting a diet.

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**CH:** MediaShop AG | Industriering 3 | 9491 Ruggell | Liechtenstein
**EU:** MediaShop GmbH | Schneiderstraße 1, Top 1 | 2620 Neunkirchen | Austria